



bobo motion 2.0

User Guide


Sep 2023



- 1 Download the **Bobo Pro 2.0** app from Google Play or the App Store.
- 2 Register within the app. Please note that the email you use for registration will be the address to which all reports will be sent.
- 3 Enter your product activation code.
You can find it printed on the side of the motion sensor.
- 4 Ensure that your Bluetooth and location settings are turned on, and grant all app requested permissions.
- 5 To power on the motion sensor, press the power button briefly.
- 6 Follow the app's guidance to properly position the Velcro straps on your patient and attach the sensor.

To replace the battery:

1. unscrew the safety screw.
2. pull the battery drawer in the arrow direction.
3. insert new battery CR2032 3V.
4. close the drawer and fasten the screw back.



What's in the box?

Motion sensor

Straps



- 30cm X1
- 60cm X1
- 90cm X1

Download app



Power On/Off

- Power on: short press**
Blue -> Orange
- Searching mode**
Blue blinking
- Connected**
Green blinking
- Power off - long press**
Wait for 2 Red blinkings

Specifications

Wireless connection to tablets and smartphones
 LED & audio indicators
 Sensors: 3-axis accelerometer, Gyroscope
 Operating temp. : -5°C to 45 °C (23 °F to 113 °F)
 Battery Type: Coin CR2032 3V - **Non-Rechargeable**
 Battery Life: Approx 400 hours of continuous use
 Communication module: Bluetooth low energy (BLE)

BO&BO Ltd / info@bobo-balance.com
 +1 888-300-3153

Manufactured in China
 E.A.R. - Obelis S.A. Boulevard Général_Wahis 53, B-1030 Brussels, BELGIUM.



"Do not ingest battery, Chemical Burn Hazard"

This product contains a coin/button cell battery. If the coin/button cell battery is swallowed, it can cause severe internal burns in just 2 hours and can lead to death.
 Keep new and used batteries away from children.
 If the battery compartment does not close securely, stop using the product and keep it away from children.
 If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.