



Open your mind to

BrainWear[®]

The most accurate, cost-effective whole brain measuring technology on the market

EMOTIV is the pioneer of simple to use, high quality, mobile EEG technologies. Recognized and validated by the scientific community and developers worldwide, our technology has been cited in over 4000 publications to date.



EMOTIV | EPOC Flex[™]

EPOC Flex is our highest density and most flexible wireless EEG system for research professionals.

Combines the award-winning wireless technology of our EPOC+ headset with the flexibility and high density afforded by more traditional EEG head cap systems.

With a subscription to our **PRO License**, the EPOC+ provides access to dense array, high quality, raw EEG data.



High density

EPOC Flex is a wireless control box that works alongside the EasyCap system. Record from any of the standard 10-20 EEG positions for up to 32 channels.



Maximum flexibility

Choose from different cap sizes to accommodate a wide range of head sizes and support awake or sleep studies with two control box positions.



Professional mobile EEG

EPOC Flex transmits data wirelessly at 128 Hz, so you can record high resolution brain data without being tethered to a computer.



Saline or gel

Flex Saline features EPOC+ felt based sensors that can be hydrated on or off the head. Traditional gel sensors offer the highest signal quality.