


# EmotivPRO

---


EmotivPRO is compatible with the following operating systems:

- Windows 10 (64-bit) v1607+
- macOS 10.13 High Sierra or above

EmotivPRO requires a minimum of 8 GB RAM to operate. For Windows, please ensure your system has a minimum specification of Core i5-3xxx processor or higher.

 For optimal rendering performance we advise to use a maximum resolution of 1920 x 1080. Downscaling the resolution on 4k, 5k displays helps to significantly improve performance.

 We recommend setting the screen scaling option to 100%.

 Some performance issues are expected if you run EmotivPRO and EmotivBCI concurrently on the same machine.