

Kardia Mobile ECG by AliveCor

Our smallest mobile ECG measurement device (8 by 3 cm, FDA clearance and CE certificate) is suitable for everyday use. See if your ECG is within the norm or if you are affected by ventricular fibrillation. For iOS and Android. Easy to use, suitable for senior citizens.

EAN: 040232360961
MSRP: 142,-€



Rest easy with reliable ECG recordings that instantly let you know how your heart is doing. The 4th generation Kardia Mobile ECG by AliveCor is FDA cleared, CE marked and is the most clinically proven mobile ECG monitor on the market. Just hold the Kardia Mobile ECG by AliveCor in your hands or on your chest for 30 seconds to measure your heart's electrical activity anytime, anywhere - no probes or wires necessary. The Kardia Mobile ECG by AliveCor will wirelessly communicate with the free Kardia ECG app to let you know if your ECG is normal or if atrial fibrillation (AF) is detected. You can easily share your ECG recordings with your doctor. Track symptoms like palpitations and shortness of breath, habits like caffeine or alcohol consumption, and activities like exercise or sleep, to help your doctor get a better understanding of what may be impacting your heart health. Works with Apple iOS and most Android smartphones and tablets. Requires download of the free Kardia ECG app, available on the App Store and Google Play.

BENEFITS

- **Tried, tested, true:** get accurate electrocardiogram (ECG) readings from your smartphone or tablet
- **FDA cleared, CE marked:** clinically proven, used by leading cardiologists
- **Unlimited storage:** free unlimited ECG recordings
- **Instant analysis:** AliveCor Instant Analysis tells you when your ECG is normal and detects atrial fibrillation (AF) in your ECG
- **Simple and easy:** one-time setup and activation
- **4th Generation:** New and improved sensors and app

THE KARDIA APP - THE HEART JOURNAL

It's never been easier to track your symptoms, activities and diet to help you and your doctor figure out what's causing your heart rhythm disturbances

- **Journaling Made Easy:** After each ECG, simply tap on a Symptom, Activity or Diet tag to add it to your recording
- **Create Entries Between ECGs:** Enter a tag or quick note to keep track of how you're feeling, what you've eaten and your activities at any time, even without recording an ECG
- **Beat Fluctuation (BFx):** With each ECG, see how your heartbeat changes from beat to beat

THE PROVIDER DASHBOARD

The Provider Dashboard is a free web-based application to help you review your patients' ECG data. Once you connect to a patient in the Provider Dashboard, their future recordings automatically start showing up. Then, you can quickly review the metrics that you deem most important (e.g. symptoms, ECG analysis reports, heart rate).

TECHNICAL INFORMATION

- Harder, more durable, scratch resistant electrodes
- Textured electrode surface, making it easier to ensure proper finger placement for recording
- Rigorously tested and meets our highest quality standards
- **Size:** (L) 3.3in x (W) x 1.4in x (D) 0.2 in,
- **Weight:** 0.6oz (includes 3.0V, CR2016 battery)

Apple Compatibility:

- iPhone 6 Plus / 6s Plus (iOS 9.2), iPhone 6 / 6s (iOS 9.2), iPhone 5 / 5s / 5c (iOS 9.2), iPhone 4s (iOS 9.2), iPod Touch 5G (iOS 9.2), iPad Air / Air 2 (iOS 9.2), iPad Mini (iOS 9.2)

Android Compatibility:

- LG Nexus 5 (Android 6.0.1), Samsung Note 5 (Android 6.0.1), Samsung Note 3 (Android 5.x), Galaxy S6 / S5 / S4 (Android 5.x), Galaxy S3 (Android 4.x), HTC One M7 (Android 4.x), Jitterbug Touch 3 / 2 (Android 4.x)
- **Warranty:** <https://www.alivecor.com/warranty>

Apps:



Kardia App

Accessories:



Kardia Mobile EKG Set



Attachment plate



Carry Pod

EAN: 4260497571018

EAN: 040232262852

EAN: 850214007087

The Kardia products match the European Directive 93/42/CEE in accordance with the German Medical Products Law and the Medical Products Regulation. They are classified as medical products (CE 0459) accordingly and fulfill all requirements regarding quality and manufacturing.

The Kardia Mobile has a CE approval and is available for purchase in the following countries: USA, Australia, Bahrain, Bangladesh, Belgium, Canada, Chile, France, Germany, Hong Kong, Hungary, Ireland, India, Israel, Italy, Jamaica, Kuwait, Netherland, Luxembourg, Malta, Norway, Pakistan, Poland, Spain, Switzerland, Trinidad and Tobago, New Zealand, Oman, Qatar, United Arab Emirates and U.K.

In all other countries it is exclusively for medical Professional staff.