



## Diagnostic ECG at your fingertips.


Kardia Mobile produces a diagnostic quality ECG on a smart phone or tablet in just **30 seconds**.

It's quick and easy to use with no gels, wires or cables and it detects atrial fibrillation, a common arrhythmia that can be the underlying cause of stroke.

The recorded ECG's can be reviewed, emailed, saved and printed.

Kardia Mobile is a reliable and accurate device that's recommended and used by clinicians around the world (Class IIa medical device and FDA cleared).

Offering the ability to quickly rule out abnormalities in patients with palpitations, Kardia saves time and costs, accelerating diagnosis and therapy.




Kardia mobile



Place two fingers from each hand on each of the sensors to produce an ECG of lead I. A live ECG trace is shown on screen and it takes **just 30 seconds** to make a recording and produce a rhythm strip.

Works with both iOS and Android phones and tablets.

ECG recording does not require Wifi or Bluetooth.



Instant analysis provides accurate detection of atrial fibrillation. Cardiologists reported AF and atrial flutter **sensitivity of 100%** in the ECG rhythm strip.<sup>1</sup>

Monitor patients post ablation and post CABG. A fast and economical tool for screening patients. Detect problems early by screening routinely for AF – reduce time to therapy.

Check patients with palpitations to rule out arrhythmias.

Use as part of a routine health check.

Take on home visits for a quick and simple ECG check.



Patients can capture events **whenever or wherever** they happen. Kardia Mobile need never leave their side when attached to a smart phone. When symptoms occur, a patient can instantly record the event and send it to their clinician by email.

Free download of the Kardia app from iTunes or Google Play.

Additional features available with Premium service subscription (30 day free trial included).

1. Tarakji KG, Wazni OM, Callahan T, Kanj M, Hakim AH, Wolski K, et al. Heart Rhythm. 2015;12(3):554-9. NICE MIB35 after the full stop.



**ORDER NOW AT: [Alivecor.com](http://Alivecor.com) or via our customer support line**

AliveCor UK Tel: +44 (0) 333 301 0433 email [uksupport@alivecor.com](mailto:uksupport@alivecor.com)

AliveCor and Kardia Mobile are trademarks and trade names of AliveCor, Inc. in the United States.

v. KMBUK276172