



## MShirt 3 – Technical Specification

2020-09-17

<b>Name</b>	MShirt 3 (part of Mbody product family)
<b>Type</b>	#43-0383 MB ... #43-0883 MB
<b>Used with</b>	MCell 3
<b>For applications</b>	Muscle Monitor 3 (for Windows) ErgoLink (for Windows)
<b>Channels</b>	8 EMG biosignal channels (4 on left side, 4 on right side):  Wrist extensors & flexors (combined) Biceps Deltoid group Trapezius
<b>Textile materials</b>	Carvico Morea® 80% Polyamide micro, 20% Lycra (170g/m2) Carvico Spider® 80% Polyamide micro, 20% Lycra (150g/m2) Electrodes silver coated yarns
<b>Washing</b>	Fully machine washable at 30°C (86°F)   <p>Avoid exposing sensors to solvent.</p>
<b>Operating temperature</b>	0°C ... +40°C (32°F ... + 104°F)
<b>Humidity</b>	10 ... 90 % RH (non-condensing)
<b>Water resistance</b>	Protected against water splashes during use. Waterproof during washing. Let dry carefully before use with MCell 3.
<b>Sizes</b>	S, M, L, XL, XXL (Unisex)
<b>Weight</b>	~ 430 g (15,2 oz)
<b>Recycling</b>	 <p>This product contains electronic components and materials. Product must be disposed as electrical and electronic equipment.</p>
<b>Standards and approvals</b>	Oeko-Tex® Standard 100 CE, FCC, CMIIT, RoHS, WEEE