

QardioBase

QardioBase is the first step to a happier, healthier you. It combines smart features, ease of use and elegant design to make weight and health management an integral part of your everyday life.



EAN: 854894005173

MSRP: 179,-€

QardioBase measures more than just your weight. It tracks your body mass index (BMI) and analyzes your full body composition including muscle mass, bone and water composition and body fat. Set goals and reminders, track your progress with charts and graphs. Choose your preferred mode to reach your ideal weight by turning your phone into your personal coach. QardioBase automatically recognizes each family member so everybody can have a truly unique experience.

With its minimalist curved design, this smart scale will be an elegant addition to your home. QardioBase is sleek enough to blend into any room and beautiful enough to be the center of attention. With its discreet display, QardioBase was made for design-conscious living, complementing the most elegant of environments. It provides precise measurement on any surface, be it a marble bathroom floor or a cozy soft carpet.

BENEFITS

- Weight & Full Body Composition
- Smart Feedback
- Pregnancy Mode
- Multiple Users

FUNCTION

1. **Weight & More:** QardioBase smart scale measures your weight, BMI (body mass index), and full body composition (body fat %, muscle, water and bone mass).
2. **Pregnancy Mode:** QardioBase tracks your progress during this happy time and beyond, measuring your progress week-by-week during and after pregnancy.
3. **Made for sharing:** QardioBase body fat scale is designed for the whole family to enjoy, cleverly recognizing multiple users of the same household automatically.
4. **Minimalist elegance:** QardioBase body composition scale is sleek

enough to blend into any room, beautiful enough to be the center of attention and accurate on any surface.

5. **Smart Feedback:** QardioBase offers this unique feedback mode, using smile, rather than numbers to track your progress towards your individual goal.
6. **Haptic Feedback:** With haptic feedback, a gentle vibration lets you know your measurement has been taken, suitable for those with visual impairment.
7. **Powerful Data sharing:** Share your measurement data with family and friends automatically or send your readings to your doctor via email, using the Qardio App.
8. **Qardio connected Ecosystem:** QardioBase Wi-Fi scale integrates with the Qardio ecosystem using Qardio App - allowing users to store and share their data via Qardio's secure cloud.

TECHNICAL INFORMATION

- **Measures & Calculates:** Weight, BMI, % body fat, muscle bone, water composition
- **Modus:** Normal, Weight Only (safe mode for those with pace makers) Smart Feedback, Pregnancy
- **Weight:** 4 lbs (1.8kg) including eight AAA batteries
- **Dimension:** 1.95 cm Hx 32 cm W
- **Power Source:** 8 x 1.5V batteries (Size AAA), supplied. Lasts 1 year with normal use
- **Materials:** Glass top with ABS bottom housing
- **Measurement Range:** 9-396 Lbs (5-180 Kg)
- **Units:** Kg, Lbs, St
- **Wireless Sync:** WI-FI(802.11 (b/g/n) compatible with WEP/ WPA/WPA2 security protocols, Bluetooth 4.0
- **Requires:** Bluetooth 4.0
- **Measurement accuracy:** 0.2 Lbs (100g) with 4 Sensors
- **Display LED:** Displays Weight, % body fat, user selection, expressive icons, Haptic Feedback
- **Works with:** Apple iOS7.0 or later, Apple Watch